

SHORE DUROMETER

- Tablet SHORE durometer for different scales (A, C, D, F, 00)
- Functions: Go/NoGo, Max/Min, Formula, Timer, Linear correction, Resolution selecting
- Memory manager: 2000 values, Folders system, Statistics mode, Memory data transfer
- **4 MODES DATA TRANSFER:** WIRELESS to MDS app (Windows, Android, iOS, MacOS); WIRELESS HID, WIRELESS HID+MAC, USB HID
- Calibration certificate included (ISO17025 (Ilac MRA))



Item No	Shore Scale	Range	Resol	Indenter			Force	Preset	Go/NoGo	Max/Min	Formula	Timer	Temp comp	Linear corr	Calibr date	FW update	Recharge	Memory	WIRELESS	USB
				Diameter	Indenting foot	Extension														
		unit	unit	mm		mm	N													
21100101	A	0-100	0,01	1,40	35° truncated cone D=0,79mm	2,54	8,05	•	•	•	•	•	•	•	•	•	•	Folders / 2000 val	Long range / HID	USB HID / Update
21100102	C				30° cone R=0,1mm		44,45	•	•	•	•	•	•	•	•	•				
21100103	D			1,20 mm spherical R	1,111		•	•	•	•	•	•	•	•	•					
21100104	00			2,40	0,55-4.3		•	•	•	•	•	•	•	•	•					
21100105	F																			

MULTISCALE SHORE DUROMETER




INTERCHANGEABLE INDENTERS

MADE IN UKRAINE



Item No	Shore Scale	Indenters pcs
21100110	A,C	2
21100111	A,D	2
21100120	A,C,D	3
21100100	A,C,D,00	4
21100130	A,C,D,00, F	5

OPTIONAL ACCESSORIES



IoT MDS CONNECT DISPLAY UNIT
USB, WI-FI, RJ-45, RS-485, LoRa OUTPUT